

PRO-Clean™

Rapid Protein Residue Test

PRO-Clean is a quick and easy way to accurately monitor the cleanliness of surfaces to help ensure product quality. PRO-Clean detects protein residues left on a surface after cleaning. Simply swab a surface, release the reagent and if protein residue is present the reagent will turn purple. The color change provides a semiquantitative measure of the surface cleanliness. The more contamination present, the quicker the color change to purple and the darker the color. PRO-Clean quickly validates the hygiene of a surface, allowing immediate corrective action to be taken if necessary.



DETECTS PROTEIN RESIDUES

Key Features

- Results in under 10 minutes
- All-in-one device — no instrumentation
- 18-month shelf life at refrigerated temperatures (2 - 25°C)
- Easy to interpret
- Affordable

NO INSTRUMENTATION


Ideal for

- Meat processors
- Food processors
- Restaurants and food service
- Grocery stores
- Industrial applications

RESULTS IN

LESS THAN 10 MIN

Premoistened Swab



Americas

941 Avenida Acaso
Camarillo, CA 93012
USA

Tel: 805.388.8007

Fax: 805.388.5531

Email: enquiries@hygienea.net

Web: www.hygienea.net

Europe

Unit 11 Wenta Business Centre
Colne Way Wafford
Hertfordshire WD2 4ND UK

Tel: +44 (0)1923 818821

Fax: +44 (0)1923 818825

Email: enquiries@hygienea.net

Web: www.hygienea.net

ORDER # PRO-100

(100 tests per box)