



CERTIFICATE OF ACCREDITATION

ANSI-ASQ National Accreditation Board

500 Montgomery Street, Suite 625, Alexandria, VA 22314, 877-344-3044

This is to certify that

Northland Laboratories – Green Bay
1030 Parkview Road
Green Bay, WI 54304

has been assessed by ANAB
and meets the requirements of international standard

ISO/IEC 17025:2005

while demonstrating technical competence in the field of

TESTING

Refer to the accompanying Scope of Accreditation for information regarding the types of tests to which this accreditation applies.

AT-1490
Certificate Number


ANAB Approval

Certificate Valid: 12/19/2017-10/20/2018
Version No. 003 Issued: 12/19/2017



This laboratory is accredited in accordance with the recognized International Standard ISO/IEC 17025:2005. This accreditation demonstrates technical competence for a defined scope and the operation of a laboratory quality management system (refer to joint ISO-ILAC-IAF Communiqué dated April 2017).



SCOPE OF ACCREDITATION TO ISO/IEC 17025:2005

Northland Laboratories – Green Bay

1030 Parkview Road
 Green Bay, WI 54304
 Krysta Renier
 920-336-7465
 krenier@northlandlabs.com

TESTING

Valid to: **October 20, 2018**

Certificate Number: **AT-1490**

Microbiological

Specific Tests and/or Properties Measured	Specification, Standard, Method, or Test Technique	Items, Materials or Product Tested	Key Equipment or Technology
Aerobic Plate Count	FDA/BAM Ch 3	Food	Plate
Aerobic Plate Count	CMMEF Ch 3 & 8	Environmental	Plate
Aerobic Plate Count by Petrifilm	AOAC 990.12	Food	Petrifilm
Aerobic Plate Count by Petrifilm	CMMEF Ch 3 & 8	Environmental	Petrifilm
Coliform/ <i>E. coli</i> Count by VRBA	FDA/BAM Ch 4	Food	Plate
Coliform/ <i>E. coli</i> Count by VRBA	CMMEF Ch 3 & 9	Environmental	Plate
Coliform/ <i>E. coli</i> Count by MPN	FDA/BAM Ch 4	Food	Multiple Tube
Coliform/ <i>E. coli</i> Count by MPN	CMMEF Ch 3 & 9	Environmental	Multiple Tube
Coliform/ <i>E. coli</i> Count by Petrifilm	AOAC 991.14 / 998.08	Food	Petrifilm
Coliform/ <i>E. coli</i> Count by Petrifilm	CMMEF Ch 3 & 9	Environmental	Petrifilm
<i>Enterobacteriaceae</i> by Petrifilm	CMMEF Ch 3 & 9	Food / Environmental	Petrifilm
<i>Staphylococcus aureus</i> Count, coagulase positive	FDA/BAM Ch 12	Food	Plate
<i>Staphylococcus aureus</i> Count, coagulase positive	CMMEF Ch 3 & 39	Environmental	Plate



Microbiological

Specific Tests and/or Properties Measured	Specification, Standard, Method, or Test Technique	Items, Materials or Product Tested	Key Equipment or Technology
<i>Staphylococcus aureus</i> Count, coagulase positive by Petrifilm	AOAC 2003.07	Food	Petrifilm
<i>Staphylococcus aureus</i> Count, coagulase positive by Petrifilm	CMMEF Ch 3 & 39	Environmental	Petrifilm
Yeast/Mold Count	FDA/BAM Ch 18	Food	Plate
Yeast/Mold Count	CMMEF Ch 3 & 21	Environmental	Plate
Yeast/Mold Count by Petrifilm	AOAC 997.02	Food	Petrifilm
Yeast/Mold Count by Petrifilm	CMMEF Ch 3 & 21	Environmental	Petrifilm
Rapid Yeast/Mold Count by Petrifilm	3M Petrifilm 6475/6477 (AOAC PTM 121301)	Food	Petrifilm
Rapid Yeast/Mold Count by Petrifilm	AFNOR № 3M 01/13-07/14	Environmental	Petrifilm
Listeria genus by ELFA	AOAC 2004.06	Food	BioMerieux VIDAS
Listeria genus by ELFA	Listeria spp, VIDAS 30700 (AOAC RI 981202)	Environmental	BioMerieux VIDAS
Listeria genus by PCR	iQ Check Listeria spp. 730003459 (AFNOR № BRD 07/13-08/07)	Food / Environmental	Bio-Rad iQCheck
<i>Listeria monocytogenes</i> by ELFA	AOAC 2004.02	Food	BioMerieux VIDAS
<i>Listeria monocytogenes</i> by ELFA	AFNOR № BIO 12/9-07/02	Environmental	BioMerieux VIDAS
<i>Listeria monocytogenes</i> by PCR	iQ Check Listeria monocytogenes 730003592 (AFNOR № BRD 07/10-04/05)	Food / Environmental	Bio-Rad iQCheck
Salmonella genus by ELFA	AOAC 2011.03	Food / Environmental	BioMerieux VIDAS
Salmonella genus by PCR	iQ Check Salmonella 730003613 (AFNOR № BRD 07/06-07/04)	Food / Environmental	Bio-Rad iQCheck

Microbiological

Specific Tests and/or Properties Measured	Specification, Standard, Method, or Test Technique	Items, Materials or Product Tested	Key Equipment or Technology
<i>E. coli</i> O157:H7 by PCR	iQ Check E.coli 0157:H7 730003460 (AFNOR № BRD 07/15-06/08)	Raw Beef, 375g	Bio-Rad iQCheck
Residual Aerobic Plate Count	CMMEF Ch 3	Packaging Film	Plate
Residual Coliform Plate Count	CMMEF Ch 3	Packaging Film	Plate

Chemical

Specific Tests and/or Properties Measured	Specification, Standard, Method, or Test Technique	Items, Materials or Product Tested	Key Equipment or Technology
Fat (Babcock)	SMEDP 15.082	Dairy	Babcock Bottles
Salt	SMEDP 15.053	Dairy	Corning Chloride Analyzer
pH	AOAC 981.12	Dairy	pH meter
Moisture	SMEDP 15.111	Dairy	Vacuum oven
Aflatoxin (B1, B2, G1, G2)	Neogen Veratox #8030 (AOAC RI 050901)	Food	Stat Fax 4700
Gluten (Gliadin R5) Allergen	Neogen Veratox #8510 (AOAC RI 061201)	Food	Stat Fax 4700
Milk Allergen	Neogen Veratox #8470	Food	Stat Fax 4700
Soy Allergen	Neogen Veratox #8410	Food	Stat Fax 4700

Note:

1. This scope is formatted as part of a single document including Certificate of Accreditation No. AT-1490.



Vice President